

Periprosthetic Fxs:

Hip

Utku KANDEMIR, MD

Professor

Dept of Orthopaedic Surgery
University of California, San Francisco






Questions

- How to know if the stem loose or not?
- How to Fix?
- MIPO vs ORIF?
- When to use an allograft strut?
- Are cables necessary?

Vancouver classification of hip periprosthetic fractures

Vancouver classification relies on:

1. The level of the fracture
2. If the prosthesis is stable or not
3. the quality of the bone



Hip periprosthetic fractures	
Type A	Peritrochanteric fractures
	AG: greater trochanter
	AL: lesser trochanter
Type B	Around or just below the tip of the femoral stem
	B1: stable stem
	B2: loose stem
	B3: loose implant with substantial bone loss
Type C	fractures occur well below the implant



Treatment Goal

- Return to previous function
- Early Mobilization



- ***REDUCTION:*** Anatomic Alignment
- Fracture Union
- Stable Prosthesis



LT fx , GT fx









GT fx - Osteolysis



GT fx - Osteolysis
















NEXT GENERATION
CABLE-CLAW
LOCKING SCREW
PLATES



GT Fx

- Plate: Screw + Claw+ Cable




Critical Q: Stem *Stable* vs *Loose*



Hip periprosthetic fractures	
Type A	Peritrochanteric fractures
	AG: greater trochanter
	AL: lesser trochanter
Type B	Around or just below the tip of the femoral stem
	B1: stable stem
	B2: loose stem
	B3: loose implant with substantial bone loss
Type C	fractures occur well below the implant


Critical Q: Stem *Stable* vs *Loose*



Stem Stability



Stem Loose: *Revision*



B1: Stable Implant: ORIF



Stem: Stable vs Loose

How to check?

- Preinjury imaging – if available
 - Subsidence, haloes
 - X-rays and CT scan
 - Implant bone interface proximally
-
- Intraoperative
 - Through Fracture
 - Through ARTHROTOMY



Stem: Stable vs Loose: *Implications*

- ORIF vs Revision Arthroplasty + Fixation
- Planning:
 - Implants
 - Surgeon(s)
 - Positioning
 - *Plan B, Plan C*

Fixation / ORIF:



What you think?



Cable only fixation proximally






What you think?



- ❑ INADEQUATE proximal fixation
- ❑ SHORT construct




Around the Stem


UNICORTICAL LOCKING SCREWS AND
CABLES +/- bicortical/transcortical
screws

- *Screws:* Keep length and anti-torsion
- *Cables:* anti pull out/anti-bending stress





84 yo





Tangential Bicortical Locked Fixation Improves Stability in Vancouver B1 Periprosthetic Femur Fractures: A Biomechanical Study


Gregory S. Lewis, PhD, Cyrus T. Caroom, MD, Hwabok Wee, PhD, Darin Jurgensmeier, MD, Shane D. Rothermel, MD, Michelle A. Bramer, MD, and John Spence Reid, MD






Tangential Bicortical Locked Fixation Improves Stability in Vancouver B1 Periprosthetic Femur Fractures: A Biomechanical Study


Gregory S. Lewis, PhD, Cyrus T. Caroom, MD, Hwabok Wee, PhD, Darin Jurgensmeier, MD, Shane D. Rothermel, MD, Michelle A. Bramer, MD, and John Spence Reid, MD

- Cable only constructs half as strong or less





6 months
f/u




Fixation / ORIF Strategies

- ❑ REDUCE it well and BUILD it strong
- ❑ *SCREWS ONLY* are stronger than *CABLES ONLY*
- ❑ SPAN the whole bone



ORIF: OPEN vs MIPO




How to Fix

- Plate
 - Plate + Allograft strut
 - Plate + Plate
- +**
- CABLES around the STEM




Needed Secondary Implant







Case example

89 year old female

- Ground level fall
- Remote history of Right THA





88 yo







88 yo, 2 months postop

□ WBAT



7 months Postop




Revision




1 month postop Revision 8 months post-index Sx




Revision #2



6 weeks postop Revision#2




4 months postop Revision#2




10 months postop Revision#2



Revision #3: 24 months post-index



2 yr post revision #3



Secondary Fixation

- Allograft Strut

- Plate




Secondary Fixation: Strut vs Plate

▣ **STRUT:**

- Increase bone stock
- Add to stability of fixation +
- Technically *harder*

▣ **PLATE:**

- -Add to stability to fixation +++
- Technically *easier*





Allograft Strut






a



1 yr f/u





When to use STRUT?

- ❑ Bone DEFECT / Bone loss
- ❑ OSTEOPOROTIC bone w very thin cortices
- ❑ As a secondary implant, NOT isolated!!!!




CABLES

- At least 2 cables around the stem
- Metallic > Synthetic > Suture
- Better around the whole fixation (plates or plate+strut)



Fixation: How Long?


- ❑ Bone quality
- ❑ Fracture pattern



- ❑ Better to span the whole Femur

65 yo

65 yo, cemented THA





65 yo, cemented THA





69 yo

69 yo





69 yo




69 yo, 1yr f/u




83yo


83 yo



83 yo



83 yo



62 yo



62 yo





62 yo



62
yo




62 yo, 6 months f/u




57 yo


57 yo




57 yo



57 yo



57 yo, f/u



Take Home Messages

- Around Stem: UNICORTICAL LOCKING SCREWS AND CABLES +/- bicortical/transcortical screws
- REDUCE it well and BUILD it strong
- Consider 2 implants unless strong Bone
- SPAN the whole bone



THANK YOU



Utku.Kandemir@ucsf.edu

