

Removing Well Fixed Implants in THA

UCSF/ Hoag Arthroplasty for the Modern Surgeon Course 2025


Erik Hansen, MD

Professor of Clinical Orthopaedic Surgery

Division Chief, Adult Reconstruction

William R. Murray Endowed Chair in Orthopaedics

Sept 27, 2025



Disclosures

- CORIN- Consultant
- Smith & Nephew- Advisory Board

OUTLINE

- Principles
- Femoral Implant Removal
- Acetabular Component Removal
- Summary



PRINCIPLES

- Preoperative planning is critical
 - Prior op reports, templating (planned ETO), instrument/implant request
 - Cup removal- OD/ID, liners for explant; femoral heads for retained stem
 - Stem removal- proprietary extraction tool; liners for retained cup
- Minimize iatrogenic bone loss
 - Maximize success of reconstruction
- Patience, patience, patience
 - Yet know when you need to switch gears (e.g. ETO)

FEMUR



WORKFLOW

- Expose lateral shoulder of implant
- Expose proximal femur-implant interface circumferentially
- Flexible osteotomes
- Extraction device

Femoral Implants

- **Single wedge taper stems**
- Modular primary stems
- Dual wedged/ Triple tapered collared stem
- Fully HA compaction broach type stem
- Cemented stems
- Diaphyseal engaging stems



Femoral Implants

- Single wedge taper stems
- **Modular primary stems**
- Dual wedged/ Triple tapered collared stem
- Fully HA compaction broach type stem
- Cemented stems
- Diaphyseal engaging stems



Femoral Implants

- Single wedge taper stems
- Modular primary stems
- **Dual wedged/ Triple tapered collared stem**
- Fully HA compaction broach type stem
- Cemented stems
- Diaphyseal engaging stems



Femoral Implants

- Single wedge taper stems
- Modular primary stems
- Dual wedged/ Triple tapered collared stem
- **Fully HA compaction broach type stem**
- Cemented stems
- Diaphyseal engaging stems



Femoral Implants

- Single wedge taper stems
- Modular primary stems
- Dual wedged/ Triple tapered collared stem
- Fully HA compaction broach type stem
- **Cemented stems**
- Diaphyseal engaging stems




Femoral Implants

- Single wedge taper stems
- Modular primary stems
- Dual wedged/ Triple tapered collared stem
- Fully HA compaction broach type stem
- Cemented stems
- **Diaphyseal engaging stems**



“Clearing the Shoulder”

- Do This In All Cases!
 - Especially varus remodeling, troch overgrowth, and long straight stems
- Allows access to lateral interface
 - Subsequent burring/ osteotomes
 - Extraction handle threads
- Facilitates revision stem reaming trajectory





Exposing the Proximal Bone Implant Interface

- Pencil tipped burr circumferentially
 - Vs. double sided recip saw for A-P sides
- Creates safe bone-implant interface
- Helps avoid the osteotomes from skiving and creating perforations
- Tips
 - Be wary of lateral shoulder and angle
 - Irrigate & suction continuously
 - Don't hesitate to remove some proximal bone



Flexible Osteotomes



- A-P: Straight blades
- Laterally: Curved blade
- Medially: Reverse curve



Flexible Osteotomes

- Tip-

- Measure blade vs length of stem fixation surface
- Ensure bevel side of blade adjacent to implant, and strike plate ergonomically placed for retrograde mallet blows



Universal Femoral Extraction Set




SHUKLA H/P

Universal Hip Implant Extraction Solution

and
MOD SHUKLA H/P
Universal Hip Implant Extraction Solution

Systems 5 & 6 of 15



STRIKE PLATE FRAME


e applies crucial extraction force in-line with the Harmonics help break up bone ingrowth, making helping reduce the risk of fractures.



The wide strike plate allows surgeons to hit with a large degree of accuracy, reducing the chance of mishits.



Taking out Modular Neck Stems



Gotta Collar?? No Problem


- Burr collar w/ metal cutting tip
- Resect some calcar to give yourself space



Exodus


Revision Hip System
for Femoral Stem Removal

Surgical Technique



WATSON EXTRACTION SYSTEM™


A Novel Surgical Instrument Set for the Extraction of Well-Fixed Cementless Stems in Revision Total Hip Arthroplasty



Curved lateral blades (CLAT)



Straight lateral blades (SLAT)




Straight medial blades (SMED)






Fully HA Coated Compaction Broaching Stem

- Additionally described techniques
 - Smooth K wires under fluoro
 - Hitting distal before proximal
- Tip
 - Be careful of perforations
 - Low threshold for ETO






Cemented Stems

- Stems are smooth so relatively easy to tamp out
 - May need to just remove the most proximal cement
- Left with a cement mantle to remove
 - Moreland set +/- burr
 - Work circumferentially w/ radial splits
- Tip
 - Scrutinize your XRs to understand cement mantle
 - Clear proximal before going distal



Distal Cement Plug/Mantle


- Strategy depends on quantity and quality of distal cement mantle
 - Drill & Tap vs back hoe
 - Ultrasonic device (Oscar, Ultradrive)
 - ETO
- Tip
 - Imperative to get all cement out before reaming for rev stem



Distal Cement Plug/Mantle



- Strategy depends on quantity and quality of distal cement mantle
 - Drill & Tap vs back hoe
 - Ultrasonic device (Oscar, Ultradrive)
 - ETO
- Tip

 - Imperative to get all cement out before reaming for revision stem



Extended Trochanteric Osteotomy- Indications

- **Removal of well-fixed femoral components* /cement**
- Improving acetabular exposure
- Facilitating implantation of revision components
- Correcting deformity (varus remodeling)
- Tip-
 - *Episiotomy may suffice



Extended Trochanteric Osteotomy- Planning



- Long enough to successfully get components out

- Short enough to keep an intact tube

- Allow for at least 4cm of scratch fit distal to ETO


- Tip-
 - Make sure you get full length femur films to template

- Measure twice (ruler & fluoro) and cut once



Extended Trochanteric Osteotomy

- Distal limb rounded with burr




Extended Trochanteric Osteotomy

- Start the back side with an osteotome




Extended Trochanteric Osteotomy

- *Lever open with multiple wide osteotomes*



Extended Trochanteric Osteotomy

- The finished product...




Extended Trochanteric Osteotomy

- Gigli saw medially



Cylindrical Diaphyseal Engaging Stems

- Cut stem at proximal aspect of cylindrical stem
- Trephine distal fragment



Conical Splined Revision Stem

- ETO
- Cannot use trephines because of tapered distal geometry
- Helpful instruments
 - Pencil tip burr
 - Curved flexible osteotomes
 - Smooth Steinman pins for splines
- Tip
 - Remove proximal body if modular
 - Have proprietary extraction device for distal stem



One stage Exchange

ACETABULUM






WORKFLOW

- Expose the acetabular component circumferentially
- Remove liner
- Remove screws
- Explant system

PE Liner Removal



- Good ole' osteotome at rim
- Fancy drill/ screw
- Tip-
 - Know idiosyncratic liner locking mechanism
 - Try not to destroy liner as can be used for explant




Zimmer Trilogy

Metal or Ceramic Liner Removal

- Proprietary extraction device
- Bone tamp on rim of cup
 - Creates vibration
- Suction cup



ZB



Smith Nephew





Figure 18



Stryker



Figure 47
Alternative Bearing Extractor





Figure 48
Placement of Alternative Bearing Extractor

Depuy

Removing Screws


- Know how many you need to go after
- Know screw head
 - Hex- Most systems
 - Star- Stryker Trident 1 pre-2018
 - Cruciform- Depuy SROM cup
- Stripped screw head?
 - Metal cutting burr the head
- Broken screw shaft?
 - Take proximal out through cup; trephine remainder later





Broken Screw Removal Set



Exposing Bone Implant Interface





Acetabular Gouges




Old school technique- sometimes helpful for non-hemispheric cups, and those w/ spikes

Explant Cup Out System


- Need to know cup OD/ID
 - If you don't have prior Op Report...
 - Preop: Set ratio of actual/measured head/cup on XR
 - Intra-op: Many cups have info inscribed on rim
- Short/ Long blades
- Work the cup systematically
 - Think “clock face”, start at 12 o'clock
- Tip-
 - Have scrub tech organized ahead of time
 - Careful around the femoral retractor
 - Evaluate your preop Xray for challenges





Large Monoblock MOM Cups

- Standard Explant system heads only go to 36mm
- Though you can use gouges
- Better yet- Use either trial bipolar heads or DM heads



Automatic Acetabular Extractors



Depuy Kincise2



Brasseler EZX



Stryker EZout



SUMMARY

- Removing well fixed components can be exhausting and time consuming
 - Preop planning is critical
 - Numerous instruments now available
 - Nothing takes place of proper technique
-